

SKC特殊採樣介質儲存方式及保存期限一覽表

| 儲存方式     | 需低溫冷藏(-8°C以下)之產品             |  | 不需低溫冷藏(但應放置蔭涼處)之產品                |   |
|----------|------------------------------|--|-----------------------------------|---|
| 保存期限     |                              |  |                                   |   |
| 一個月      | 225-9014                     | 預裝二片玻璃纖維濾紙(nitrite-impregnated)之濾匣                             |                                   |   |
|          | 226-165A                     | 矽膠(mercuric cyanide)吸附管  |                                   |   |
| 三個月      | 225-9010                     | 預裝二片玻璃纖維濾紙(veratrylamineredyl& di-n-ocyl phthalate)之濾匣         |                                   |   |
|          | 225-9013                     | 預裝一片玻璃纖維濾紙與纖維素墊片(1-(2-pyridyl)piperazine)之濾匣                   |                                   |   |
|          | 225-9018                     | 預裝一片玻璃纖維濾紙與一片聚脂濾紙(mercuric chloride)之濾匣                        |                                   |   |
|          | 225-9019                     | 預裝二片玻璃纖維濾紙(2,4-dinitrophenyl-hydrazine and phosphoric acid)之濾匣 |                                   |   |
|          | 225-9020                     | 預裝三片玻璃纖維濾紙(2,4-dinitrophenyl-hydrazine and phosphoric acid)之濾匣 |                                   |   |
| 225-9021 | 預裝二片玻璃纖維濾紙(veratrylamine)之濾匣 |  |                                   |   |
| 六個月      | 225-9002                     | 預裝一片玻璃纖維濾紙與纖維素墊片(1-(2-pyridyl)piperazine)之濾匣                   |                                   |   |
|          | 225-9009                     | 預裝二片玻璃纖維濾紙(1-(2-pyridyl)piperazine)之濾匣                         |                                   |   |
|          | 225-9011                     | 預裝一片25mmAE玻璃纖維濾紙(methoxyphenyl piperazine)之濾匣                  |                                   |   |
|          | 500-100                      | UME-100 被動式採樣器(Formaldehyde 採樣用)                               |                                   |   |
| 八個月      | 225-9007                     | 預裝一片玻璃纖維濾紙(mercuric acetate)之濾匣                                |                                   |   |
| 一年       | 225-9003                     | 預裝二片玻璃纖維濾紙(2,4-dinitrophenyl-hydrazine and phosphoric acid)之濾匣 | 225-9005                          | 預裝MCE 預濾片與墊片; 纖維素濾紙與墊片(Sodium carbonate)之濾匣               |
|          | 225-9015                     | 預裝二片超厚玻璃纖維濾紙(ascorbic acid)之濾匣                                 | 225-9006                          | 預裝Teflon 預濾片與多孔塑膠墊片; 特殊潔淨銀膜濾紙與多孔塑膠墊片之濾匣                   |
|          | 226-27                       | XAD-2 (2-Hydroxymethyl piperidine)吸附管                          | 225-9016                          | 預裝二片PVC與墊片(Sodium formate) / Metrical filter(uncoated)之濾匣 |
|          | 226-30-07                    | XAD-2 (p-Anisidine)吸附管   |                                   |   |
|          | 226-38-03                    | 石油基活性炭(Hydrobromic acid)吸附管                                    |                                   |   |
|          | 226-54                       | XAD-2 (2-Hydroxymethyl piperidine)吸附管                          |                                   |   |
|          | 226-68                       | JXC 活性炭 Drierite(Hydroquinone)吸附管                              |                                   |   |
|          | 226-70A                      | 矽膠 (p-Methoxyphenol)吸附管  |                                   |   |
|          | 226-117                      | XAD-2 (2-Hydroxymethyl piperidine) 吸附管                         |                                   |   |
|          | 226-118                      | XAD-2 (2-Hydroxymethyl piperidine)吸附管                          |                                   |   |
|          | 226-151                      | 活性炭 (Proprietary coating)吸附管                                   |                                   |   |
|          | 226-153                      | XAD-2 (di-n-butylamine)吸附管                                     |                                   |   |
|          | 500-400                      | UME-400 被動式採樣器(Amines 採樣用)                                     |                                   |   |
|          | 575-005                      | 被動式採樣器(Treated Anasorb PC) (環氧乙烷採樣用)                           |                                   |   |
| 二年       | 226-40                       | 分子篩(Triethanolamine)/Oxidizer吸附管(2 tubes);                     | 225-9001                          | 預裝玻璃纖維濾紙(1-(2-pyridyl)piperazine)之濾匣                      |
|          | 226-119                      | 矽膠(2,4-Dinitrophenylhydrazine) 吸附管                             | 225-9004                          | 預裝玻璃纖維濾紙(Sulfuric acid)之濾匣                                |
|          | 226-119A                     | 矽膠(2,4-Dinitrophenylhydrazine)吸附管                              | 225-9012                          | 預裝玻璃纖維濾紙(Sulfuric acid)之濾匣                                |
|          | 226-120                      | 矽膠(2,4-Dinitrophenylhydrazine) 吸附管 with built-in ozone scruber | 225-9022                          | 預裝一片PTFE濾紙與一片玻璃纖維濾紙(MAMA)之濾匣                              |
|          |                              |  | 226-10-06                         | 矽膠 吸附管(Sulfuric acid)                                     |
|          |                              |  | 226-23                            | XAD-2 (Octanoic acid) 吸附管                                 |
|          |                              |  | 226-29                            | Anasorb 747 (Sulfuric acid) 吸附管                           |
|          |                              |  | 226-30-16                         | XAD-2/Glass Fiber Filter(OVS) 吸附管                         |
|          |                              |  | 226-30-16A                        | XAD-2 / Glass Fiber Filter (OVS) 吸附管                      |
|          |                              |  | 226-30-18                         | XAD-2 (Naphthylisothiocyanate) 吸附管                        |
|          |                              |  | 226-35-01                         | Tenax 吸附管   |
|          |                              |  | 226-35-02                         | Tenax /Tenax 吸附管  |
|          |                              |  | 226-40-02                         | Molecular Sieve(Triethanolamine) 吸附管                      |
|          |                              |  | 226-42                            | Silica Gel (Sulfuric acid) 吸附管                            |
|          |                              |  | 226-42-02                         | Firebrick (Gas chrom-R) (Sulfuric acid) 吸附管               |
|          |                              |  | 226-53                            | Silica Gel (Sulfuric acid) 吸附管                            |
|          |                              |  | 226-55                            | Silica Gel (Sodium hydroxide) 吸附管                         |
|          |                              | 226-56   | Tenax/Glass Fiber Filter(OVS) 吸附管 |   |

SKC特殊採樣介質儲存方式及保存期限一覽表

| 儲存方式 | 不需低溫冷藏(但應放置蔭涼處)之產品 |  |            |                               |
|------|--------------------|--|------------|-------------------------------|
| 保存期限 |                    |  |            |                               |
| 二年   | 226-57             | XAD-7/Glass Fiber Filter (OVS) 吸附管             | 226-339    | Tenax TA 吸附管                  |
|      | 226-57A            | XAD-7/Glass Fiber Filter (OVS) 吸附管             | 226-340    | Tenax TA 吸附管                  |
|      | 226-58             | XAD-2/Quartz Filter (OVS) 吸附管                  | 226-341    | Anasorb CMS 吸附管               |
|      | 226-58A            | XAD-2/Quartz Filter (OVS) 吸附管                  | 226-345    | TenaxGR/Anasorb GCB1 吸附管      |
|      | 226-61             | Silica Gel / Charcoal (Sodium hydroxide) 吸附管   | 226-346    | Anasorb GCB1/ Anasorb CMS 吸附管 |
|      | 226-63             | Silica Gel 吸附管                                 | 226-347    | Anasorb GCB2/GCB1/CMS 吸附管     |
|      | 226-67             | Anasorb CSC (Potasium hydroxide) 吸附管           | 226-348    | TenaxGR/Anasorb GCB1 吸附管      |
|      | 226-73             | Anasorb CSC (t-Butylcatechol) 吸附管              | 226-349    | Anasorb GCB1/ Anasorb CMS 吸附管 |
|      | 226-80             | Anasorb 747 (Potasium hydroxide) 吸附管           | 226-350    | Anasorb GCB2/GCB1/CMS 吸附管     |
|      | 226-92             | PUF 吸附管(低流量用)                                  | 226-356    | Anasorb GCB1 吸附管              |
|      | 226-96             | XAD-7 ([NBD] Chloride) 吸附管                     | 226-357    | Tenax TA 吸附管                  |
|      | 226-98             | XAD-7 (Phosphoric acid) 吸附管                    | 226-358    | Chromosorb 106 吸附管            |
|      | 226-99             | Silica Gel/Glass Fiber Filter(OVS) 吸附管         | 226-359    | Anasorb GCB1 吸附管              |
|      | 226-104            | Chromosorb 102 吸附管                             | 226-360    | Tenax TA 吸附管                  |
|      | 226-109A           | Chromosorb 104 吸附管                             | 226-361    | Chromosorb 106 吸附管            |
|      | 226-116A           | Porapak-T / Porapak-T 吸附管                      | 520-02A    | 被動式採樣器吸附劑夾囊(汞蒸氣採樣用)           |
|      | 226-124            | PUF / Tenax / PUF 吸附管(低流量用)                    | 575-001    | 被動式採樣器(Charcoal)(有機蒸氣採樣用)     |
|      | 226-125            | Glass beads (2,4,6-Trichlorophenol)/Silica 吸附管 | 575-001A   | 被動式採樣器(Charcoal)(有機蒸氣採樣用)     |
|      | 226-126            | PUF / Glass Fiber Filter 吸附管(低流量用)             | 575-001B   | 被動式採樣器(Charcoal)(有機蒸氣採樣用)     |
|      | 226-129            | PUF / XAD-2 / PUF 吸附管(高流量用)                    | 575-001C   | 被動式採樣器(Charcoal)(有機蒸氣採樣用)     |
|      | 226-131            | PUF 吸附管(高流量用)                                  | 575-001MC  | 被動式採樣器(Charcoal)(二氯甲烷採樣用)     |
|      | 226-133            | Anasorb747 吸附管                                 | 575-001MCA | 被動式採樣器(Charcoal)(二氯甲烷採樣用)     |
|      | 226-134            | Tenax TA 吸附管                                   | 575-002    | 被動式採樣器(Anasorb 747)(有機蒸氣採樣用)  |
|      | 226-135            | Tenax TA / Petroleum Charcoal 吸附管              | 575-002A   | 被動式採樣器(Anasorb 747)(有機蒸氣採樣用)  |
|      | 226-142            | Carbon Beads(KOH) /Teflon Filter 吸附管           | 575-002B   | 被動式採樣器(Anasorb 747)(有機蒸氣採樣用)  |
|      | 226-143            | PUF / XAD-2 / PUF 吸附管(低流量用)                    | 575-002C   | 被動式採樣器(Anasorb 747)(有機蒸氣採樣用)  |
|      | 226-302            | Anasorb CMS 吸附管                                | 575-003    | 被動式採樣器(Anasorb 727)(有機蒸氣採樣用)  |
|      | 226-323            | Tenax TA 吸附管                                   | P22692R    | 再充填PUF 吸附管(低流量用)              |
|      | 226-330            | Anasorb CMS / GCB1 / GCB2 吸附管                  | P226129R   | 再充填PUF/XAD-2/PUF 吸附管(高流量用)    |
|      | 226-336            | Tenax TA 吸附管                                   | P226131R   | 再充填PUF 吸附管(高流量用)              |
|      | 226-337            | Anasorb CMS 吸附管                                | P226143R   | 再充填PUF/XAD-2/PUF 吸附管(低流量用)    |
|      | 226-338            | Anasorb CMS / GCB1 / GCB2 吸附管                  |            |                               |